

Your diabetes data, simply there.




With the mySugr® app it's quick and easy to collect all your important diabetes therapy data in one place through a growing number of connected devices, integrations and manual logging.

mySugr is created by people with diabetes for people with diabetes and eases the complexity of your daily diabetes routine. We envision a world where a person with diabetes can live a full and healthy life, freed from worrying about the daily grinds of their therapy. With over 4 million registered app users in more than 70 countries, we are well on our way to achieving our goal to make diabetes suck less®!

mySugr Report

Download and access all your diabetes data in just a few clicks! The mySugr Report compiles your data for easy printing or digital sharing. Thanks to the Report function, you can access your data anytime, anywhere, and even send an overview directly to your doctor and healthcare team.



Better check-ups with better info!



- Easily share your diabetes information with multiple reporting formats¹
- Easy to spot patterns: View your report directly from your smartphone for quick decision-making
- Share your report electronically with your healthcare provider between appointments
- Custom reports showcase just the data you want to share by selecting your time frame of choice for PDF- & Excel Report.

¹ mySugr Reports are available in CSV-format with mySugr Basic for downloading all your diabetes data or in PDF and Excel format with mySugr PRO, incl. the option to choose a time frame of your choice (2 weeks, 3 months or Custom for your PDF Report).



Take telehealth to the next level!




Why wait three months to update therapy settings? Send an overview of your diabetes data directly to your healthcare team any time. Help them make informed, constructive decisions about your unique diabetes therapy and put yourself back in the driver's seat! mySugr Reports provide you with:

- **Blood glucose (BG) average for the week**
- **Standard deviation for the week:** how much your blood sugar has gone above or below your average glucose – the “swing” factor)
- **Activity time for the week**
- **The percentage of hypos:** the number of blood sugar entries below your target range
- **The percentage of hypes:** the number of blood sugar entries above your target range



- 1 Open the mySugr app and select “Reports” from the menu
- 2 Create your personal mySugr Report in CSV format with mySugr Basic or in PDF or Excel with mySugr PRO
- 3 Select your file format, date range, and other options. The mySugr Report in PDF or Excel format is detailing 2 weeks, 3 months or a custom date range of your diabetes data
- 4 Select “Export”
- 5 Choose “Share” (print, email, message, or save)²



² When you share your report (e.g. via email), your personal data leaves the mySugr app. We cannot guarantee the privacy of your data in this regard.


Connect your Accu-Chek® blood glucose meter and auto-activate mySugr PDF Reports for free



Get started and let mySugr help you get the diabetes results you need!

Connect your Accu-Chek blood glucose meter and automatically import your blood sugar results. Did you know? When connecting your meter, you unlock mySugr PRO for free to get advanced app features like PDF Reports and so much more! Here is how it works:

- 1 Download the mySugr app
- 2 Unlock mySugr PRO by connecting your Accu-Chek meter
- 3 Create your personal mySugr Report PDF or Excel with mySugr PRO
- 4 Share your PDF Report



Frequently featured in popular diabetes media, (Medical Tribune, The Diabetes Times, The Guardian and Deutsche Ärztezeitung) mySugr is the app that simplifies daily diabetes therapy. mySugr eases the process of logging and features a home screen where you can view graphs, daily stats, data, and much more! With mySugr, a host of features come together in one easy to use tool.

