ACCU-CHEK® Guide

EN

Control Solutions

Cat. No. 07748906001

Control Level 1

Control Level 2

For use with the Accu-Chek Guide, the Accu-Chek Guide Link, and the Accu-Chek Guide Me blood glucose monitoring systems.

For in vitro diagnostic use only.

Before You Start Testing

Purpose

The control solution acts like blood when you apply it to an Accu-Chek Guide test strip. Performing a control test lets you know the meter and test strips are working properly.

You should perform a control test when:

- you open a new test strip box.
- you left the test strip container open.
- you think the test strips are damaged.
- you want to check the meter and test strips.
- the test strips were stored in extreme temperatures, humidity, or both.
- you dropped the meter.
- your test result does not match how you feel.
- you want to check if you are performing the test correctly.

WARNING

- Do not put the control solutions in your mouth. Do not inject the control solutions. They are for use outside the body only.
- Choking hazard. Small parts. Keep away from children under the age of 3 years.

About the control solutions

- Use only the Accu-Chek Guide control solutions.
- If the bottle is open or damaged before using the control solution for the first time, if the cap is not fully closed, if you see any damage to the cap or bottle, or if anything prevents the cap from closing properly, do not use the control solution. Contact the Accu-Chek Customer Care Service Center at 1-800-858-8072.
- Close the control solution bottle tightly after use.
- The unopened control solution is stable until the **Use By** date printed on the bottle label (24 months from date of manufacture) if the control solution is stored properly. Do not use control solution that is past the **Use By** date.
- Write the date you open the bottle on the bottle label. The control solution must be discarded 3 months from that date or on the **Use By** date printed on the bottle label, whichever comes first.
- Store the control solution bottle at room temperature (39–86 °F) and 10–90 % relative humidity. Do not freeze the control solution.
- The meter automatically recognizes the difference between the control solution and blood.
- The control solution can stain fabric. Remove stains by washing with soap and water.

STEP 1 Getting Ready to Test

You need the Accu-Chek Guide, the Accu-Chek Guide Link, or the Accu-Chek Guide Me meter, the Accu-Chek Guide test strips and a Level 1 or Level 2 bottle of the Accu-Chek Guide control solution.

NOTE

- The control level is printed on the bottle label.
- The acceptable ranges for the control solutions are printed on the test strip container.
- Check that the control solutions and test strips are not expired.

STEP 2 Performing a Control Test


1. Insert the metallic end of the test strip into the meter. The meter turns on. Place the meter on a flat surface, like a table.
2. Check the **Use By** date on the test strip container.

Do not use test strips past the **Use By** date.



3. Select the control solution you want to test, either Level 1 or Level 2. You will enter the level later in the test.
4. Remove the control bottle cap. Wipe the tip of the bottle with a tissue. **On the Accu-Chek Guide and the Accu-Chek Guide Link meters:**

5. When **Apply drop** appears on the display, squeeze the bottle until a tiny drop forms at the tip of the bottle. Touch the drop to the yellow edge of the test strip. Do not put control solution on top of the test strip. **Analyzing** appears when there is enough control solution in the test strip. Wipe the tip of the bottle with a tissue, then cap the bottle tightly.

6. **Control Result** and the control bottle symbol appear. Do not remove the test strip yet. Press  or  to select the control level you tested.

7. Press . **Within range** and  appear if the control result is within range. **Out of range** and  appear if the control result is out of range.

On the Accu-Chek Guide Me meter:

5. When you see the flashing  appear on the display, squeeze the bottle until a tiny drop forms at the tip of the bottle. Touch the drop to the yellow edge of the test strip. Do not put control solution on top of the test strip. The flashing  appears when there is enough control solution in the test strip. Wipe the tip of the bottle with a tissue, then cap the bottle tightly.

6. The control result, the control bottle symbol, and the flashing **L1** or **L2** appear on the display. Press  to confirm the control level you tested.

Note: Press  or  to alternate between **L1** and **L2**.

7. **OK** appears if the control result is within range. **Err** appears if the control result is out of range.

You can also compare your control result to the range printed on the test strip container.

STEP 3 Understanding Control Results

In Range Control Results

The meter and test strips are working properly if **Within range** and  or **OK** appear on the display, or if the control result is within the range printed on the test strip container.

Out of Range Control Results

Check the table below to see what may be causing the problem if **Out of range** and  or **Err** appear on the display, or if the control result is outside of the range printed on the test strip container.

Troubleshooting Checks	Action
Were the test strips or control solutions expired?	Discard the test strips or control solution if either is past the Use By date. If the control solution was opened more than 3 months ago, discard it. Repeat the control test with an unexpired test strip and an unexpired control solution.
Did you wipe the tip of the control solution bottle before use?	Wipe the tip of the bottle with a tissue. Repeat the control test with a new test strip and a fresh drop of control solution.
Were the caps on the test strip container and the control solution bottle always closed tightly?	Replace the test strips or control solution if you think either was uncapped for some time. Repeat the control test.
Was the test strip used immediately after it was removed from the test strip container?	Repeat the control test with a new test strip and a fresh drop of control solution.
Were the test strips and control solutions stored in a cool, dry place?	Repeat the control test with a properly stored test strip or control solution.
Did you follow the directions?	Review the testing steps. Repeat the control test with a new test strip and a fresh drop of control solution.
Did you choose the correct control solution level, either 1 or 2, when you performed the test?	If you chose the wrong control solution level, you can still compare the control result to the range printed on the test strip container.
Are you still unsure of the problem?	Contact the Accu-Chek Customer Care Service Center at 1-800-858-8072.

Technical Information

Each bottle contains glucose, buffer, salts, non-reactive ingredients, preservative, FD & C Blue #1.

Explanation of Symbols

	Follow instructions for use
	Global Trade Item Number
	Manufacturer

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© 2024 Roche Diabetes Care
ACCU-CHEK, ACCU-CHEK GUIDE, and ACCU-CHEK GUIDE ME are trademarks of Roche.
Roche Diabetes Care, Inc.