



July 10, 2017

To Whom It May Concern:


**Subject: SELF DECLARATION**


**Blood Glucose Meter ACCU-CHEK®Guide Supplied by Roche Diabetes Care.**

Herewith we certify that the above meter for the purpose of monitoring glucose concentration in human blood are tested according to Environmental Conditions and Test Procedures for Airborne Equipment (RTCA/DO-160E) Section 21 Emission of Radio Frequency Energy.

When the meter is in flight mode, the meter fulfills the requirements according to this standard, category M.

The meter does not affect the performance of aircrafts.

  
Karen Roberts  
Director of Verification and Validation  
Roche Diabetes Care Inc.

  
John Krug  
Director of Quality  
Roche Diabetes Care Inc.